

The Multi-Task Genie™ Family

THE ORIGINAL



Scientific Industries' reputation for high quality and reliability is world renowned - to date, more than 1,000,000 genuine Vortex-Genie mixers have been used in any type of laboratory! We continue our tradition as the "Gold Standard" in vortexing and mixing with our complete Multi-Task Genie family of vortex mixers.

Vortex-Genie® 1

- Violent single high-speed vortexing
- Easy-to-use Touch On operation
- Large, 6mm orbit for aggressive vortexing of difficult samples
- May be used with Pop-Off Cup (shown) or optional 3-inch platform
- Vortex-Genie quality at an even more economical price!

Vortex-Genie® 2

- OUR MOST POPULAR MODEL
- Variable speed for precise mixing from gentle to vigorous
- Broad range of attachments for most tubes, plates and other containers
- Hands-Free or Touch On control
- May be used in cold rooms and incubators

Vortex-Genie® 2T

- Same great features as the Vortex-Genie 2
- Timer function:
 - Touch On: 1-60 sec.
 - Hands-Free: 1-60 min. or continuous
- Unit automatically stops at the end of user-programmed time
- May be used in cold rooms and incubators

U.S. Pats. 4,781,487 / 5,707,861

Ordering Information

Catalog No.	Description
SI-0136	Vortex-Genie 1, 120V, 60Hz
SI-0146	Vortex-Genie 1, 230V, 50Hz - No Plug
SI-0156	Vortex-Genie 1, 230V, 50Hz - Euro Plug
SI-0166	Vortex-Genie 1, 230V, 50Hz - British Plug
SI-0176	Vortex-Genie 1, 230V, 50Hz - Swiss Plug
SI-0186	Vortex-Genie 1, 100V, 50/60Hz
SI-0236 (Model G560)	Vortex-Genie 2, 120V, 60Hz
SI-0246 (Model G560E)	Vortex-Genie 2, 230V, 50Hz - No Plug
SI-0256 (Model G560E)	Vortex-Genie 2, 230V, 50Hz - Euro Plug
SI-0266 (Model G560E)	Vortex-Genie 2, 230V, 50Hz - British Plug
SI-0276 (Model G560E)	Vortex-Genie 2, 230V, 50Hz - Swiss Plug
SI-0286	Vortex-Genie 2, 100V, 50/60Hz
SI-T236	Vortex-Genie 2T, 120V, 60Hz
SI-T246	Vortex-Genie 2T, 230V, 50Hz - No Plug
SI-T256	Vortex-Genie 2T, 230V, 50Hz - Euro Plug
SI-T266	Vortex-Genie 2T, 230V, 50Hz - British Plug
SI-T276	Vortex-Genie 2T, 230V, 50Hz - Swiss Plug
SI-T286	Vortex-Genie 2T, 100V, 50/60Hz

Specifications-All Vortex-Genie mixers

Weight	4 Kg (8.8 lb)
Base Dimensions (D x W x H)	165 x 122 x 165mm 6.5 x 4.8 x 6.5in
Supplied	- Vortex-Genie 1 - Vortex-Genie 2* - Vortex-Genie 2T*
Timer	- Vortex-Genie 2T Touch On: 1-60 seconds Hands-Free: 1-60 minutes or continuous

* See back page for additional Vortex-Genie 2 and Vortex-Genie 2T accessories.



ETL Listed Conforms to UL Std. 3101-1
ETL Listed Certified to CAN/CSA C22.2



Scientific Industries, Inc.
Home of the Multi-Task Genie™ Family
70 Orville Drive, Bohemia, NY 11716 U.S.A.
Phone: (631) 567-4700 • Fax: (631) 567-5896
Toll Free: 888-850-6208
Customer Service: custsvc@scientificindustries.com
www.scientificindustries.com

Multi-Task Genie™ Family Accessories

Our Vast Range Of Accessories Allows For The Ultimate In Multi-Task Versatility For Any Application!




Fig. 1
Model H301
Multiple Sample
Starter Set




Fig. 4
Model H302
Large Sample Set




Fig.	Catalog No.	Description	Fig.	Catalog No.	Description
1	504-0233-00	Model H301 Multiple Sample Starter Set	4	503-0282-00	Model H302 Large Sample Set
--	146-6005-00	6-inch Platform (not shown separately)	5	504-0039-00	Recessed Platform
2	504-0234-00	Microtube Foam Insert (2)	6	503-0278-00	9-16mm Tube Foam Inserts (2)
3	504-0235-00	Microplate Foam Insert (2)	7	503-0279-00	14-29mm Tube Foam Inserts (2)
11	146-3011-00	Pop-off Cup	8	503-0280-00	29-37mm Tube Foam Inserts (2)
12	OK-0500-902	3-inch Platform with Rubber Cover	9	503-0277-00	Blank Foam Inserts (2)
13	SI-0563	TurboMix™ Attachment, 1.5ml Tubes	--	503-0283-00	Assorted Foam Inserts (4)
--	SI-0562	TurboMix Attachment, 2.0ml Tubes (not shown)	10	568-0001-00	Elastic Bands (2)
--	0A-0563-010	1.5ml Microtube Holder* (not shown separately)	--	SI-0510	Two-Tier Microplate Foam Insert (not shown)
14	0A-0563-011	2.0ml Screw-cap Microtube Holder*	15	SI-0511	Large Ampule/Tube Attachment
			16	SI-0570	Small Ampule/Tube Attachment

* May be used with TurboMix attachment or Disruptor Genie

NOTE: Vortex-Genie 1 may only be used with pop-off cup and 3-inch platform

U.S. Pat. 5,707,861



Scientific Industries, Inc.
Home of the Multi-Task Genie™ Family
70 Orville Drive, Bohemia, NY 11716 U.S.A.
Phone: (631) 567-4700 • Fax: (631) 567-5896
Toll Free: 888-850-6208
Customer Service: custsvc@scientificindustries.com
www.scientificindustries.com