

ROTO-SHAKE GENIE®

Rock, Rotate, or 3-D Orbit Any Container!



More multi-task products from the makers of the Vortex-Genie®!

- Combines rocking and rotating in one compact and versatile unit
- Rocking tray accommodates large open containers
- Magnetic strips secure sealed bags for completely uniform soaking and washing
- Optional 3-D orbital attachment transforms the unit into an orbital shaker
- Easily cleaned or decontaminated
- **MAGNETIC** stainless steel platform may be easily set to rock or rotate almost any container
- A variety of attachments are available to provide maximum application versatility
- Small footprint and economical price complement any type of laboratory
- May be used in cold rooms and incubators between 0-38°C

U.S. Pat. 5,918,979

Ordering Information

Catalog No.	Description
SI-1100	Roto-Shake Genie, 120V
SI-1101	Roto-Shake Genie, 230V, No Plug
SI-1102	Roto-Shake Genie, 230V, European Plug
SI-1103	Roto-Shake Genie, 230V, British Plug
SI-1104	Roto-Shake Genie, 230V, Swiss Plug
SI-1106	Roto-Shake Genie, 240V, Australian Plug
SI-1105	Roto-Shake Genie, 100V

Specifications

Supplied	1 Clip Plate each for 10-13 mm, 15-17mm, 28-30mm Tubes, 4 Bag Mounting Strips, Non-Skid Tray
Rotating	5-35 RPM
Rocking	10-70 cycles per minute
Rocking Angle	±10 degrees from horizontal
Platform Capacity	4.5 Kg (10 lb)
Weight	11.8 Kg (26 lb)
Dimensions (D x W x H)	252 x 508 x 267mm (9.9 x 20 x 10.5in)
Depth with Rocker Tray	305mm (12in)
Environmental	Installation Category 2, 0° - 38°C (32° - 100°F), 95% max. relative humidity




Home of the Multi-Task Genie™ Family
 70 Orville Drive, Bohemia, NY 11716 U.S.A.
 Phone: (631) 567-4700 • Fax: (631) 567-5896
 Toll Free: 888-850-6208
 Customer Service: custsvc@scientificindustries.com
scientificindustries.com

Roto-Shake Genie[®] Accessories

A Vast Range Of Accessories For The Ultimate Multi-Task Versatility In Any Application!



Accessories attach magnetically to platform!



Accessories supplied with the Roto-Shake Genie

Fig.	Cat. No.	Description	Fig.	Cat. No.	Description
1	SI-1123	1 Clip Plate for 12 each 10-13mm tubes	8	SI-1127	Expansion Kit - two stainless steel trays, 24 bag mounting strips
2	SI-1124	1 Clip Plate for 6 each 15-17mm tubes	9	SI-1130	2 Clip Plates for Hybridization Tubes
3	SI-1125	1 Clip Plate for 3 each 28-30mm tubes	10	SI-1134	2 Universal Clip Plates
4	SI-1120	6 Clip Plates for 12 each 10-13mm tubes	11	SI-1135	1 Magnetic Covered 100 Microtube Box
5	SI-1121	6 Clip Plates for 6 each 15-17mm tubes	12	SI-1136	1 Clip Plate for 100ml Volumetric Flasks
6	SI-1122	6 Clip Plates for 3 each 28-30mm tubes	13	0A-1100-070	Non-Skid Tray
7	SI-1126	Bag Mounting Strips (pack of 16)	14	SI-1150	3-D Orbital Shaker Attachment

Plasticware, glassware and other consumable products shown are not included.



Home of the Multi-Task Genie™ Family
 70 Orville Drive, Bohemia, NY 11716 U.S.A.
 Phone: (631) 567-4700 • Fax: (631) 567-5896
 Toll Free: 888-850-6208
 Customer Service: custsvc@scientificindustries.com
scientificindustries.com