


BAG ROTATOR GENIE™

Powerful and rugged rotator for effective mixing of transfer packs or random sealed bags, maintaining platelets and other sample components in suspension without risk of damage.



Economic Rotator For Transfer Packs and Bags

- Wire basket holds up to 40 standard 300ml transfer packs (5 packs per compartment)
- Specialized basket design prevents bag movement
- Can effectively mix contents of almost any sealed bag type
- Gentle yet thorough mixing eliminates the risk of sample damage or trauma
- Bags rotate “end-over-end” with a fixed speed of 3 RPM regardless of the number of bags
- Maintains platelets in even suspension throughout the bag
- Durable and robust motor provides reliability and reproducibility
- Heavy and durable base prevents unit from tipping over

Ordering Information

Catalog No.	Description
SI-1000	Bag Rotator Genie, 120V
SI-1001	Bag Rotator Genie, 230V, No Plug
SI-1002	Bag Rotator Genie, 230V, European Plug
SI-1003	Bag Rotator Genie, 230V, British Plug
SI-1004	Bag Rotator Genie, 240V, Swiss Plug
SI-1006	Bag Rotator Genie, 120V, Australian Plug
SI-1005	Bag Rotator Genie, 100V

Specifications

Weight	7.7 kg (17.0 lb) with basket
Base Dimensions (D x W x H)	252 x 508 x 267mm (9.9 x 20.0 x 10.5in)
Speed	Fixed rotation, 3 RPM ±0.5 RPM
Capacity	Up to (40) 300ml transfer packs or random bags
Warranty	2-Years



Home of the Multi-Task Genie™ Family
80 Orville Drive, Suite 102
Bohemia, NY 11716 U.S.A.
Phone: (631) 567-4700 • Fax: (631) 567-5896
Toll Free: 888-850-6208
Customer Service: custsvc@scientificindustries.com
scientificindustries.com