

The Multi-Task Genie™ Family



HIGH SPEED MIXER



THE ORIGINAL



INTEGRATED TIMER



DIGITAL CONTROL

Vortex-Genie® 1	Vortex-Genie® 2	Vortex-Genie® 2T	Digital Vortex-Genie® 2
<ul style="list-style-type: none"> Violent single high-speed vortexing Easy-to-use Touch On operation Large, 6mm orbit for aggressive vortexing of difficult samples May be used with Pop-off Cup (shown) or optional 3-inch Platform May be used in cold rooms and incubators 	<ul style="list-style-type: none"> OUR MOST POPULAR MODEL Variable speed for precise mixing from gentle to vigorous Broad range of attachments for most tubes, plates and other containers Hands-Free or Touch On control May be used in cold rooms and incubators 	<ul style="list-style-type: none"> Same great features as the Vortex-Genie 2 with an integrated timer Timer function: <ul style="list-style-type: none"> - Touch On: 1-60 sec. - Hands-Free: 1-60 min. or continuous Unit automatically stops at the end of user-programmed time May be used in cold rooms and incubators 	<ul style="list-style-type: none"> Same great features as the Vortex-Genie 2 with digital control and display of speed and time Digital display provides accuracy, reproducibility and repeatability Timer function: <ul style="list-style-type: none"> - Touch On: 1-99 sec. - Hands-Free: 1-99 min. or continuous May be used in cold rooms and incubators Specially configured unit for Applied Biosystems protocol available

SPECIFICATIONS

Vortex-Genie Model	Vortex-Genie 1	Vortex-Genie 2	Vortex-Genie 2T	Digital Vortex-Genie 2
Supplied	Pop-off Cup	3-inch Platform, Pop-off Cup	3-inch Platform, Pop-off Cup	3-inch Platform, Pop-off Cup
Speed (RPM)	3200 (2700 for 50Hz models)	600-3200 (2700 for 50Hz models)	600-3200 (2700 for 50Hz models)	500-3000 (2850 for 50Hz models)
Timer - Touch On: Hands-Free:	N/A	N/A	1-60 seconds 1-60 minutes or continuous	1-99 seconds 1-99 minutes or continuous
Weight	4 Kg (8.8 lb)	4 Kg (8.8 lb)	4 Kg (8.8 lb)	4 Kg (8.8 lb)
Dimensions (base) (D x W x H)	165 x 122 x 165mm (6.5 x 4.8 x 6.5in)	165 x 122 x 165mm (6.5 x 4.8 x 6.5in)	165 x 122 x 165mm (6.5 x 4.8 x 6.5in)	165 x 122 x 165mm (6.5 x 4.8 x 6.5in)
Markings / Warranty	ETL, CE / 2 years	ETL, CE / 2 years	ETL, CE / 2 years	CE / 2 years

CATALOG NUMBERS

					AB MODELS	
120V (North America)	SI-0136	SI-0236 (Model G560)	SI-T236	SI-A236	SI-A536	
230V - No Plug	SI-0146	SI-0246 (Model G560E)	SI-T246	SI-A246	SI-A546	
230V - European Plug	SI-0156	SI-0256 (Model G560E)	SI-T256	SI-A256	SI-A556	
230V - British Plug	SI-0166	SI-0266 (Model G560E)	SI-T266	SI-A266	SI-A566	
230V - Swiss Plug	SI-0176	SI-0276 (Model G560E)	SI-T276	SI-A276	SI-A576	
240V - Australian Plug	SI-0196	SI-0297 (Model G560E)	SI-T296	SI-A296	SI-A596	
100V	SI-0186	SI-0286	SI-T286	SI-A286	SI-A586	

See back page for additional Vortex-Genie 2, Vortex-Genie 2T and NEW Digital Vortex-Genie 2 accessories.

U.S. Pats. 4,781,487 / 5,707,861



Home of the Multi-Task Genie™ Family
 80 Orville Drive, Suite 102
 Bohemia, NY 11716 U.S.A.
 Phone: (631) 567-4700 • Fax: (631) 567-5896
 Toll Free: 888-850-6208
 Customer Service: custsvc@scientificindustries.com
scientificindustries.com

Multi-Task Genie™ Family Accessories

Accessories for all Vortex-Genie® 2 mixers and Vortex-Genie® Pulse



Fig. 1
Model H301
Multiple
Sample
Starter Set




Fig. 4
Model H302
Large Sample
Set



Fig.	Cat. No.	Description	Fig.	Cat. No.	Description
1	504-0233-00	Model H301 Multiple Sample Starter Set	13	SI-0564	TurboMix™ Attachment
--	146-6005-00	6-inch Platform	14	SI-0565	Universal Microtube Holder*
2	504-0234-00	Microtube Foam Insert (2)**	15	SI-0511	Large Ampule/Tube Attachment
3	504-0235-00	Microplate Foam Insert (2)**	16	SI-0570	Small Ampule/Tube Attachment
3A	504-0232-00	Deepwell Microplate Foam Insert (2)**	17	SI-H506	Horizontal 50ml Tube Holder (6 tubes)
4	503-0282-00	Model H302 Large Sample Set	18	SI-H512	Horizontal 15ml Tube Holder (12 tubes)
5	504-0039-00	Recessed Platform	19	SI-H524	Horizontal Microtube Holder (24 tubes)
6	503-0278-00	9-16mm Tube Foam Inserts (2)***	20	SI-V203	Vertical High Speed 50ml Tube Holder (3 tubes)
7	503-0279-00	14-29mm Tube Foam Inserts (2)***	21	SI-V206	Vertical High Speed 15ml Tube Holder (6 tubes)
8	503-0280-00	29-37mm Tube Foam Inserts (2)***	22	SI-V209	Vertical High Speed Microtube Holder (9 tubes)
9	503-0277-00	Blank Foam Inserts (2)***	23	SI-V506	Vertical 50ml Tube Holder (6 tubes)
--	503-0283-00	Assorted Foam Inserts (4)***	24	SI-V512	Vertical 15ml Tube Holder (12 tubes)
10	568-0001-00	Elastic Bands (2)***	25	SI-V524	Vertical Microtube Holder (24 tubes)
11	146-3011-00	Pop-off Cup	26	SI-V525	Vertical (6) 50ml, (6) 15ml, (12) Microtube Holder
12	0K-0500-902	3-inch Platform with Rubber Cover	--	0K-0236-408	Bottom Closure with Screws Kit †
--	580-2013-00	Rubber Cover for 3-inch Platform	--	SI-0510	Two-Tier Microplate Foam Insert**

* May be used with TurboMix attachment or Disruptor Genie®

** 6-inch Platform Required

*** Recessed Platform Required

† May not fit older models. Please inquire.

NOTE: Vortex-Genie 1 may only be used with Pop-off Cup and 3-inch Platform

U.S. Pat. 5,707,861



Home of the Multi-Task Genie™ Family
80 Orville Drive, Suite 102
Bohemia, NY 11716 U.S.A.
Phone: (631) 567-4700 • Fax: (631) 567-5896
Toll Free: 888-850-6208
Customer Service: custsvr@scientificindustries.com
scientificindustries.com

© 2017 Scientific Industries LTP0049 Rev. J