

The Multi-Task Genie™ Family



THE ORIGINAL



TIMER



DIGITAL

U.S. Pats. 4,781,487 / 5,707,861

Vortex-Genie Model	Vortex-Genie® 2	Vortex-Genie® 2T	Digital Vortex-Genie® 2
Supplied	3-inch Platform, Pop-off Cup	3-inch Platform, Pop-off Cup	3-inch Platform, Pop-off Cup
Speed (RPM)	600-3200	600-3200	500-3000
Timer - Touch On: Hands-Free:	N/A	1-60 seconds 1-60 minutes or continuous	1-99 seconds 1-99 minutes or continuous
Weight	4 kgs (8.8 lbs)	4 kgs (8.8 lbs)	4 kgs (8.8 lbs)
Dimensions (base) (D x W x H)	165 x 122 x 165 mm (6.5 x 4.8 x 6.5 in)	165 x 122 x 165 mm (6.5 x 4.8 x 6.5 in)	165 x 122 x 165 mm (6.5 x 4.8 x 6.5 in)
Markings / Warranty	ETL, CE / 2 years	ETL, CE / 2 years	CE / 2 years

CATALOG NUMBERS	SI-02XX	SI-T2XX	SI-A2XX	SI-A5XX
				ABI Model

Multi-Task Genie™ Family Accessories




Fig.	Cat. No.	Description	Fig.	Cat. No.	Description
2,3	504-0233-00	Model H301 Multiple Sample Starter Set	13	SI-0564	TurboMix™ Attachment
--	146-6005-00	6-inch Platform	14	SI-0565	Universal Microtube Holder - (Requires SI-0564)
2	504-0234-00	Microtube Foam Insert (2) - (Requires 146-6005-00)	15	SI-0511	Large Ampule/Tube Attachment
3	504-0235-00	Microplate Foam Insert (2) - (Requires 146-6005-00)	16	SI-0570	Small Ampule/Tube Attachment
--	504-0232-00	Deepwell Plate Insert (2) - (Requires 146-6005-00)	17	SI-H506	Horizontal 50ml Tube Holder (6 tubes)
5-10	503-0282-00	Model H302 Large Sample Set	18	SI-H512	Horizontal 15ml Tube Holder (12 tubes)
5	504-0039-00	Recessed Platform without Bands	19	SI-H524	Horizontal Microtube Holder (24 tubes)
6	503-0278-00	9-16mm Tube Foam Inserts (2) - (Requires 504-0039-00)	20	SI-V203	Vertical High Speed 50ml Tube Holder (3 tubes)
7	503-0279-00	14-29mm Tube Foam Inserts (2) - (Requires 504-0039-00)	21	SI-V206	Vertical High Speed 15ml Tube Holder (6 tubes)
8	503-0280-00	29-37mm Tube Foam Inserts (2) - (Requires 504-0039-00)	22	SI-V209	Vertical High Speed Microtube Holder (9 tubes)
9	503-0277-00	Blank Foam Inserts (2) - (Requires 504-0039-00)	23	SI-V506	Vertical 50ml Tube Holder (6 tubes)
6-9	503-0283-00	Assorted Foam Inserts (4) - (Requires 504-0039-00)	24	SI-V512	Vertical 15ml Tube Holder (12 tubes)
10	568-0001-00	Elastic Bands (2) - (Requires 504-0039-00)	25	SI-V524	Vertical Microtube Holder (24 tubes)
11	146-3011-00	Pop-off Cup	26	SI-V525	Vertical (6) 50ml, (16) 15ml, (12) Microtube Holder
12	OK-0500-902	3-inch Platform with Rubber Cover	--	OK-0236-408	Bottom Closure with Screws Kit
--	580-2013-00	Rubber Cover only for 3-inch Platform	--	SI-0513	Microplate Adaptor for Applied Biosystems