


VORTEX-GENIE® PULSE Pulsing Vortex Mixer



Vortex-Genie Pulse Accessories

More Than Just Vortexing!

- ⦿ True pulsing action provides more collisions and random motion.
- ⦿ For applications requiring more aggressive action than standard vortexing.
- ⦿ Variable "RUN/STOP" time provides quick pulsing or interval mixing.
- ⦿ Fully programmable "ON" and "REST" function reduces heat buildup within the vessel, while keeping materials in suspension. Ideal for heat sensitive samples.
- ⦿ Use with all Vortex-Genie accessories.
- ⦿ Ideal for applications such as re-suspensions of lyophilized materials, creating solutions, enzymatic reactions and much more.

Specifications

Supplied	3-inch Platform, Pop-off Cup
Speed	500 - 3000 RPM (2850 for 50Hz models)
Timer - Hands Free	1-99 minutes or continuous
Pulse - ON/OFF	0.1 seconds - 99 minutes
Weight	4 kgs (8.8 lbs)
Dimensions (D x W x H)	165 x 122 x 165mm (6.5 x 4.8 x 6.5in)



Cat. No.	Description	Cat. No.	Description
SI-P2X6	Vortex-Genie Pulse	SI-0564	TurboMix™ Attachment
504-0233-00	Model H301 Multiple Sample Starter Set	SI-0565	Universal Microtube Holder - (Requires SI-0564)
146-6005-00	6-inch Platform	SI-0511	Large Ampule/Tube Attachment
504-0234-00	Microtube Foam Insert (2) - (Requires 146-6005-00)	SI-0570	Small Ampule/Tube Attachment
504-0235-00	Microplate Foam Insert (2) - (Requires 146-6005-00)	SI-H506	Horizontal 50ml Tube Holder (6 tubes)
504-0232-00	Deepwell Microplate Foam Insert (2) - (Requires 146-6005-00)	SI-H512	Horizontal 15ml Tube Holder (12 tubes)
503-0282-00	Model H302 Large Sample Set	SI-H524	Horizontal Microtube Holder (24 tubes)
504-0039-00	Recessed Platform without Bands	SI-V203	Vertical High Speed 50ml Tube Holder (3 tubes)
503-0278-00	9-16mm Tube Foam Inserts (2) - (Requires 504-0039-00)	SI-V206	Vertical High Speed 15ml Tube Holder (6 tubes)
503-0279-00	14-29mm Tube Foam Inserts (2) - (Requires 504-0039-00)	SI-V209	Vertical High Speed Microtube Holder (9 tubes)
503-0280-00	29-37mm Tube Foam Inserts (2) - (Requires 504-0039-00)	SI-V506	Vertical 50ml Tube Holder (6 tubes)
503-0277-00	Blank Foam Inserts (2) - (Requires 504-0039-00)	SI-V512	Vertical 15ml Tube Holder (12 tubes)
503-0283-00	Assorted Foam Inserts (4) - (Requires 504-0039-00)	SI-V524	Vertical Microtube Holder (24 tubes)
568-0001-00	Elastic Bands (2) - (Requires 504-0039-00)	SI-V525	Vertical (6) 50ml, (6) 15ml, (12) Microtube Holder
146-3011-00	Pop-off Cup	OK-0236-408	Bottom Closure with Screws Kit
OK-0500-902	3-inch Platform with Rubber Cover	SI-0513	Microplate Adaptor for Applied Biosystems
580-2013-00	Rubber Cover only for 3-inch Platform		