

INCUBATOR-GENIE™

It's a Hybridization Oven, an Incubator, a Rocker/Shaker, a Rotator, ALL-IN-ONE!
Great Value & Made in the USA



Accessories Supplied with Unit

The multi-function Incubator-Genie™ is designed to heat, rotate, rock and shake all in one compact unit.

- Ideal for hybridization, digestions, or incubation
- Magnetic stainless steel platform may be easily set to rock or rotate almost any container
- Containers can be easily mounted on the magnetic platform with magnetic clips or strips
- Chamber air is gently and continuously re-circulated ensuring temperature consistency and uniformity

Supplied

Specifications

1 Clip Plate each for 10-13 mm, 15-17mm, 28-30mm Tubes, 4 Bag Mounting Strips, Non-Skid Tray, Shelf

Temperature

28°C - 75°C (82°F - 167°F)
Accuracy: ± 0.2°C, Uniformity: ± 0.5°C
Resolution: 1.0°C, Stability: ± 0.1°C

Rotating

3-35 RPM

Rocking

6-70 cycles per minute

Platform Capacity

4.5 kgs (10 lbs)

Chamber Volume

24 L (0.86 ft³)

Weight

19.5 kgs (43 lbs)

Dimensions (D x W x H)

400 x 560 x 370 mm

(Overall)

15.75 x 22 x 14.5 in

Dimensions (D x W x H)

10 x 14.25 x 10.25 in

(Chamber)

10 x 14.25 x 10.25 in

Alarms

Visible and Audible for Temperature



| Fig. | Catalog No. | Description | Fig. | Catalog No. | Description |
|------|-------------|---|------|-------------|--|
| -- | SI-140X | Incubator Genie | 9 | SI-1130 | 2 Clip Plates for Hybridization Tubes |
| 1 | SI-1123 | 1 Clip Plate for 12 each 10-13mm tubes | 10 | SI-1131 | 1 Stackable Shelf |
| 2 | SI-1124 | 1 Clip Plate for 6 each 15-17mm tubes | 11 | SI-1134 | 2 Universal Clip Plates |
| 3 | SI-1125 | 1 Clip Plate for 3 each 28-30mm tubes | 12 | SI-1135 | 1 Magnetic Covered 100 Microtube Box |
| 4 | SI-1120 | 6 Clip Plates for 12 each 10-13mm tubes | 13 | SI-1136 | 1 Clip Plate for 100ml Volumetric Flasks |
| 5 | SI-1121 | 6 Clip Plates for 6 each 15-17mm tubes | 14 | SI-1250 | 3-D Orbital Shaker Attachment |
| 6 | SI-1122 | 6 Clip Plates for 3 each 28-30mm tubes | 15 | SI-1275 | Horizontal Orbital Shaker Attachment |
| 7 | SI-1126 | Bag Mounting Strips (pack of 16) | 16 | SI-1137 | Dual Port Mixing Tube (use with SI-1130) |
| 8 | SI-1127 | Bag Expansion Kit | | | |



U.S. Pat. 5,918,979