



MULTI VORTEX-GENIE®

Large Capacity Shaker



“NEXT GENERATION” HIGH CAPACITY VORTEX MIXING

- ☉ Same true vortexing as with the original Vortex-Genie® 2.
- ☉ New ADJUSTABLE clamping system accommodates a greater variety of containers, including deep well plates and up to (4) 500mL bottles.
- ☉ New clamping system allows the user to quickly change from one type of vessel to another.
- ☉ Digital control and display of speed and time for reproducibility and accuracy, with true speed display.
- ☉ Use for gentle, orbital-style shaking or more vigorous vortexing.
- ☉ Durable and heavy zinc die casting with specially designed feet prevents “walking”.
- ☉ No additional holders or accessories necessary. Racks are included.

Ordering Information

Catalog No.	Description
SI-M2X6	Multi Vortex-Genie
SI-M010	12-Pack Multi-Tube Rack
SI-M020	12-Pack Multi-Tube Rack (set of 2)

Specifications

Supplied	Two 12-Pack Multi-Tube Racks holding up to 24 tubes of any material up to 29mm dia.
Speed	200 - 1000 RPM (up to 950 RPM for 50Hz models)
Vortexing Orbit	4mm dia. (same as Vortex-Genie 2)
Timer	1-99 minutes or continuous
Dimensions (D x W x H)	310 x 270 x 195mm (12.2 x 10.6 x 7.7in)
Weight	11.3 kgs (25 lbs)

