

QUADMAG GENIE®

Four Position Magnetic Stirrers



Analog and Digital Magnetic Stirrers - Stir up to Four Containers!

Another revolutionary product from the makers of the Vortex-Genie®!

- ⦿ Choice of TWO different models- Analog or Digital.
- ⦿ Powerful "rare earth" magnets stay coupled even under the most demanding stirring conditions.
- ⦿ Rugged ABS/Polycarbonate housing and the famous Vortex-Genie motor.
- ⦿ Digital models offer a smooth ramp-up to desired speed avoiding potential spills.
- ⦿ Digital models incorporate a timer from 1 to 99 minutes.
- ⦿ May be used in cold rooms or incubators.

Ordering Information

Catalog No.	Description
SI-4236	QuadMag Genie, 120V
SI-4246	QuadMag Genie, 230V, No Plug
SI-4256	QuadMag Genie, 230V, Euro Plug
SI-4266	QuadMag Genie, 230V, British Plug
SI-4276	QuadMag Genie, 230V, Swiss Plug
SI-4296	QuadMag Genie, 240V, Australian Plug
SI-4286	QuadMag Genie, 100V
SI-5236	Digital QuadMag Genie, 120V
SI-5246	Digital QuadMag Genie, 230V, No Plug
SI-5256	Digital QuadMag Genie, 230V, Euro Plug
SI-5266	Digital QuadMag Genie, 230V, British Plug
SI-5276	Digital QuadMag Genie, 230V, Swiss Plug
SI-5296	Digital QuadMag Genie, 240V, Australian Plug
SI-5286	Digital QuadMag Genie, 100V

Specifications

Top Surface Dimensions (D x W)	305 x 305mm (12 x 12in)
Overall Dimensions (D x W x H)	310 x 305 x 130mm (12.25 x 12 x 5.1in)
Weight	3 kgs (6.6 lbs)
Stirring Speed	250-1300 RPM
Time	1-99 min or continuous (digital models only)
Maximum Stirring Positions	4
Maximum Stirring Capacity	2 Liters Each Position (Simultaneous)
Stir Bar Maximum Length	5.1cm (2 in)
Material	ABS/Polycarbonate Non-Skid Top Pad (1)



Home of the Multi-Task Genie™ Family
 70 Orville Drive, Bohemia, NY 11716 U.S.A.
 Phone: (631) 567-4700 • Fax: (631) 567-5896
 Toll Free: 888-850-6208
 Customer Service: custsvcs@scientificindustries.com
scientificindustries.com