

ROTO-SHAKE GENIE®

ROCK, ROTATE OR 3-D ORBIT ANY CONTAINER!
Great Value & Made in the USA



- ⦿ Platform ROCKS or ROTATES
- ⦿ Magnetic Platform and clip plate system
- ⦿ Optional Stackable Tray doubles ROCKING capacity
- ⦿ May be used in cold rooms and incubators

<p>Supplied</p> <p>Rotating</p> <p>Rocking</p> <p>Platform Capacity</p> <p>Weight</p> <p>Dimensions (D x W x H) (Overall)</p> <p>Depth with Rocker Tray</p> <p>Environmental</p>	<p>Specifications</p> <p>1 Clip Plate each for 10-13 mm, 15-17mm, 28-30mm Tubes, 4 Bag Mounting Strips, Non-Skid Tray</p> <p>1-35 RPM</p> <p>2-70 cycles per minute</p> <p>4.5 kgs (10 lbs)</p> <p>11.8 kgs (26 lbs)</p> <p>252 x 508 x 267 mm</p> <p>9.9 x 20 x 10.5 in</p> <p>305 mm (12 in)</p> <p>Installation Category 2, 0° - 38°C (32° - 100°F), 95% max. relative humidity</p>
--	---



Fig.	Catalog No.	Description	Fig.	Catalog No.	Description
--	SI-110X	Roto-Shake Genie	9	SI-1130	2 Clip Plates for Hybridization Tubes
1	SI-1123	1 Clip Plate for 12 each 10-13mm tubes	10	SI-1170	1 Stackable Shelf
2	SI-1124	1 Clip Plate for 6 each 15-17mm tubes	11	SI-1134	2 Universal Clip Plates
3	SI-1125	1 Clip Plate for 3 each 28-30mm tubes	12	SI-1135	1 Magnetic Covered 100 Microtube Box
4	SI-1120	6 Clip Plates for 12 each 10-13mm tubes	13	SI-1136	1 Clip Plate for 100ml Volumetric Flasks
5	SI-1121	6 Clip Plates for 6 each 15-17mm tubes	14	SI-1150	3-D Orbital Shaker Attachment
6	SI-1122	6 Clip Plates for 3 each 28-30mm tubes	15	SI-1175	Horizontal Orbital Shaker Attachment
7	SI-1126	Bag Mounting Strips (pack of 16)	16	SI-1137	Dual Port Mixing Tube (use with SI-1130)
8	SI-1127	Bag Expansion Kit			



U.S. Pat. 5,918,979