

# SCIENTIFIC INDUSTRIES SUGGESTED CATALOG COPY

Suggested Category: Rotators

## **BAG ROTATOR—Rotate Most Bags!**

Powerful and rugged rotator for effective mixing of transfer packs or random sealed bags, maintaining platelets and other sample components in suspension without risk of damage.

- Wire basket holds up to 40 standard 300mL transfer packs (95 packs per compartment).
- Can effectively mix contents of almost any sealed bag type.
- Bags rotate “end-over-end” with a fixed speed of 3rpm regardless of the number of bags.
- Maintains platelets in even suspension throughout the bag.
- Heavy and durable base prevents unit from tipping over.

### **ORDERING INFORMATION**

<b>Catalog No.</b>	<b>Description</b>
SI-1000	Bag Rotator, 120V
SI-1001	Bag Rotator, No Plug
SI-1002	Bag Rotator, 230V European Plug
SI-1003	Bag Rotator, 230V British Plug
SI-1004	Bag Rotator, 230V Swiss Plug
SI-1005	Bag Rotator, 100V
SI-1006	Bag Rotator, 240V Australian Plug
SI-1050	Replacement Basket

### **SPECIFICATIONS**

Supplied with	Wire basket
Speed	Fixed rotation, 3 RPM +/-0.5 RPM
Dimensions	267 x 508 x 252mm (D x W x H) (10.5 x 20 x 9.9in)
Capacity	up to 40 standard 300 mL transfer packs or random bags
Weight	7.7 Kg (17 lb) with basket
Warranty	2-year manufacturer warranty