

# SCIENTIFIC INDUSTRIES SUGGESTED CATALOG COPY

Suggested Category: Mixers, Shakers, Vortexers

## **DIGITAL VORTEX-GENIE® 2 MIXER—Digital Control and Display of Speed & Time**

The Digital Vortex-Genie features the same versatility and reliability found in the original Vortex-Genie 2 along with digital control and display of speed and time

- Variable speed control allows slow shaking up to high vortexing. In addition, each attachment is designed to produce different mixing actions, allowing users to choose the correct action for each application.
- A broad range of attachments of various shapes, sizes, and materials allows for mixing of almost all common tubes or containers.
- Digital display provides accuracy, reproducibility and repeatability.
- Both Touch and Hands-Free mode automatically turns unit off after set time. Touch Mode: 1-90 seconds. Hands-Free Mode: 1-99 minutes or continuous untimed mode.
- May be used in cold rooms or incubators

### **ORDERING INFORMATION**

<b>Catalog No.</b>	<b>Description</b>
SI-A236	Digital Vortex-Genie 2, 120V
SI-A246	Digital Vortex-Genie 2, 230V, No Plug
SI-0565	Universal Microtube Holder
0K-0500-901	3-inch Platform w/o Rubber Cover
0K-0500-902	3-inch Platform with Rubber Cover
146-3011-00	Pop-off Cup
146-6005-00	6-inch Platform
503-0277-00	Blank Foam Inserts (2)
503-0278-00	9-16mm Tube Foam Inserts (2)
503-0279-00	14-29mm Tube Foam Inserts (2)
503-0280-00	29-37mm Tube Foam Inserts (2)
503-0282-00	Large Sample Set (Model H302)
503-0283-00	Assorted Foam Inserts (4)
504-0039-00	Recessed Platform
504-0233-00	Multiple Sample Starter Set (Model H301)
504-0234-00	Microtube Foam Insert (2)
504-0235-00	Microplate Foam Insert (2)
568-0001-00	Elastic Bands (2)
580-2013-00	Rubber Cover for 3-inch Platform
SI-0510	Two-Tier Microplate Foam Insert
SI-0511	Large Ampule/Tube Attachment
SI-0564	TurboMix Attachment
SI-0570	Small Ampule/Tube Attachment
SI-H506	Horizontal 50mL Tube Holder (6 tubes)
SI-H512	Horizontal 15mL Tube Holder (12 tubes)
SI-H524	Horizontal Microtube Holder (24 tubes)
SI-V506	Vertical 50mL Tube Holder (6 tubes)
SI-V512	Vertical 15mL Tube Holder (12 tubes)
SI-V524	Vertical Microtube Holder (24 tubes)
SI-V525	Vertical 50mL, 15mL, Microtube Holder
SI-V203	Vertical High Speed 50mL Tube Holder (3 tubes)
SI-V206	Vertical High Speed 15mL Tube Holder (9 tubes)
SI-V209	Vertical High Speed Microtube Holder (9 tubes)

### **SPECIFICATIONS**

Supplied with	3-inch Platform and Pop-off Cup
Dimensions (base)	165 x 122 x 165mm (D x W x H) (6.5 x 4.8 x 6.5in)
Weight	4 Kg (8.8 lb)
Warranty	2-year manufacturer warranty
Speed (RPM)	500-3000 (2850 for 50Hz models)
Timer – Touch On	1-99 seconds
Hands-Free	1-99 minutes or continuous