

SCIENTIFIC INDUSTRIES SUGGESTED CATALOG COPY

Suggested Category: Mixers, Shakers, Vortexers

VORTEX-GENIE® MAX— Heavy Duty Vortex Mixer

A fully programmable large high speed large capacity vortex mixer designed to hold and provide true vortexing in a variety of containers.

- High Speed Vortexing up to 2000 RPM.
- Holds up to 20 test tubes, vials, or cylinders of various sizes (with accessory), up to 4 microplates (with accessory) or up to six up to 250ml vessels or combination.
- Fully programmable – speed, time, reversing time, variable RUN/STOP time, and acceleration and deceleration rates.
- May be used with open containers.
- Durable and heavy zinc die casting with specially designed feet prevents “walking”.
- Ideal for repetitive protocols.

ORDERING INFORMATION

Catalog No.	Description
SI-MX200	Vortex-Genie® Max, 120V
SI-MX201	Vortex-Genie® Max, 230V, No Plug

SPECIFICATIONS

Capacity	6 250ml flasks (as supplied) 20 tubes of any kind or material up to 29mm. Dia. (50ml centrifuge tube) 4 microplates
Speed	500-2000 rpm
Timer	1 minute to 99 hours or continuous
Timer – ON/FF	.1 seconds – 99 minutes
Dimensions	310 x 270 x 195 mm (D x W x H) (12.2 x 10.6 x 7.7 in)
Weight	10.4 Kg (23 lbs)
Warranty	2-year manufacturer warranty