

SCIENTIFIC INDUSTRIES SUGGESTED CATALOG COPY

Suggested Category: Mixers, Shakers, Vortexers

NEW VORTEX-GENIE® PULSE—Pulsing Vortex Mixer For More Than Just Vortexing

A digital vortex mixer featuring true pulsing action that can be set from .1 seconds to 99 minutes for both the ON and OFF times or constant ON. Speed can also be digitally controlled. Provides true speed display.

- True pulsing action interrupts action for a more random motion and thorough mixing.
- Provides more collisions by the forming and unforming of a vortex.
- For applications requiring more aggressive action than standard vortexing but less than homogenization or aggressive disruption.
- Use with all Vortex-Genie accessories.

ORDERING INFORMATION

Catalog No.	Description
SI-P236	Vortex-Genie Pulse, 120V, 60Hz
SI-P246	Vortex-Genie Pulse, 230V, 50Hz, No Plug
SI-0565	Universal Microtube Holder for TurboMix
0K-0500-901	3-inch Platform w/o Rubber Cover
0K-0500-902	3-inch Platform with Rubber Cover
146-3011-00	Pop-off Cup
146-6005-00	6-inch Platform
503-0277-00	Blank Foam Inserts (2)
503-0278-00	9-16mm Tube Foam Inserts (2)
503-0279-00	14-29mm Tube Foam Inserts (2)
503-0280-00	29-37mm Tube Foam Inserts (2)
503-0282-00	Large Sample Set (Model H302)
503-0283-00	Assorted Foam Inserts (4)
504-0039-00	Recessed Platform
504-0233-00	Multiple Sample Starter Set (Model H301)
504-0234-00	Microtube Foam Insert (2)
504-0235-00	Microplate Foam Insert (2)
568-0001-00	Elastic Bands (2)
580-2013-00	Rubber Cover for 3-inch Platform
SI-0510	Two-Tier Microplate Foam Insert
SI-0511	Large Ampule/Tube Attachment
SI-0564	TurboMix Attachment
SI-0570	Small Ampule/Tube Attachment
SI-H506	Horizontal 50mL Tube Holder (6 tubes)
SI-H512	Horizontal 15mL Tube Holder (12 tubes)
SI-H524	Horizontal Microtube Holder (24 tubes)
SI-V506	Vertical 50mL Tube Holder (6 tubes)
SI-V512	Vertical 15mL Tube Holder (12 tubes)
SI-V524	Vertical Microtube Holder (24 tubes)
SI-V525	Vertical 50mL, 15mL, Microtube Holder
SI-V203	Vertical High Speed 50mL Tube Holder (3 tubes)
SI-V206	Vertical High Speed 15mL Tube Holder (9 tubes)
SI-V209	Vertical High Speed Microtube Holder (9 tubes)

SPECIFICATIONS

Supplied with	3-inch Platform and Pop-off Cup
Dimensions (base)	165 x 122 x 165mm (D x W x H) (6.5 x 4.8 x 6.5in)
Weight	4 Kg (8.8 lb)
Warranty	2-year manufacturer warranty
Speed (RPM)	500-3000 (2850 for 50Hz models)
Timer – Hands Free	1-99 minutes or continuous
Timer – ON/FF	.1 seconds – 99 minutes