

SCIENTIFIC INDUSTRIES SUGGESTED CATALOG COPY

NEW Genie Temp-Shaker 300™— Benchtop Orbital Shaking Incubator

Powerful Orbital Shaking up to 300RPM available in 3 different interchangeable platforms, with precise temperature control ranging from 28° Celsius to 75° Celsius. The Genie Temp-Shaker 300 allows for quick access through its unique 3-way chamber access design. Easily tilt the chamber cover left, right, or back for easy access. The super lightweight cover can be completely removed for cleaning and attaching accessories or placing containers. The unique design of temperature flow allows for minimal loss of temperature with the chamber door opened. The super-fast heaters allow the chamber to achieve a stable 75° Celsius in only 20 minutes. Never lose sight of your samples with a 360-degree view at all times.

- GTS 300 Orbital-Genie with Universal Ratcheting Clamps – provides the flexibility for using a mix of different size and shape vessels on a single platform. Unit is supplied with three ratcheting cushioned clamps. Extra ratcheting clamps for more than two rows are optional.
- GTS 300 Orbital-Genie for Standard Flask Clamps is equipped with a platform for accepting various size flasks or beakers. Clamps sold separately.
- GTS 300 Orbital-Genie with Adhering Mat provides a flat platform with adhering rubber mat perfect for use with open containers. Patented material is naturally adhering. Rinse clean to renew adhering properties.

ORDERING INFORMATION

Catalog No.

Description

SI-G1500	Genie Temp-Shaker 300 Ratcheting Clamps, 120V Standard Plug
SI-G1501	Genie Temp-Shaker 300 Ratcheting Clamps, 230V No Plug
SI-G1502	Genie Temp-Shaker 300 Ratcheting Clamps, 230V Euro Plug
SI-G1503	Genie Temp-Shaker 300 Ratcheting Clamps, 230V British Plug
SI-G1504	Genie Temp-Shaker 300 Ratcheting Clamps, 230V Swiss Plug
SI-G1505	Genie Temp-Shaker 300 Ratcheting Clamps, 100V Standard Plug
SI-G1506	Genie Temp-Shaker 300 Ratcheting Clamps, 230V Australian Plug
SI-G1600	Genie Temp-Shaker 300 Flask Clamps, 120V Standard Plug
SI-G1601	Genie Temp-Shaker 300 Flask Clamps, 230V No Plug
SI-G1602	Genie Temp-Shaker 300 Flask Clamps, 230V Euro Plug
SI-G1603	Genie Temp-Shaker 300 Flask Clamps, 230V British Plug
SI-G1604	Genie Temp-Shaker 300 Flask Clamps, 230V Swiss Plug
SI-G1605	Genie Temp-Shaker 300 Flask Clamps, 100V Standard Plug
SI-G1606	Genie Temp-Shaker 300 Flask Clamps, 230V Australian Plug
SI-G1700	Genie Temp-Shaker 300 Adhering Mat, 120V Standard Plug
SI-G1701	Genie Temp-Shaker 300 Adhering Mat, 230V No Plug
SI-G1702	Genie Temp-Shaker 300 Adhering Mat, 230V Euro Plug
SI-G1703	Genie Temp-Shaker 300 Adhering Mat, 230V British Plug
SI-G1704	Genie Temp-Shaker 300 Adhering Mat, 230V Swiss Plug
SI-G1705	Genie Temp-Shaker 300 Adhering Mat, 100V Standard Plug
SI-G1706	Genie Temp-Shaker 300 Adhering Mat, 230V Australian Plug
SI-M150	Stacking Tray
SI-1616	Double-Sided Adhesive Tape Roll, 60in. x 1in. (1.5m x25.4mm)
SI-M1617	Adhering Mat, Blue, 11.8in. x 11.8in. (300mm x 300mm)
SI-G150	GTS Non-Slip Platform
SI-G1550	GTS Ratcheting Clamp Platform
SI-G1650	GTS Flask Clamp Platform
SI-G1750	GTS Adhering Mat Platform
SI-1610	100 mL Flask Clamp (Platform accepts 16 of this size)
SI-1611	125/150 mL Flask Clamp (Platform accepts 16 of this size)
SI-1612	250 mL Flask Clamp (Platform accepts 16 of this size)
SI-1613	500 mL Flask Clamp (Platform accepts 9 of this size)
SI-1614	1 L Flask Clamp (Platform accepts 5 of this size)
SI-1615	2 L Flask Clamp (Platform accepts 3 of this size)

SPECIFICATIONS

Speed	35 – 300 RPM
Temperature Range	28 – 75 degrees Celsius
Orbit	19mm
Platform Dimensions (D x W)	305 x 305mm (12 x12in)
Overall Dimensions	450 x 480mm (17.6 x 19in)
Weight	9 Kg (20 lbs)
Chamber Height	305mm (12in)
Chamber Volume	40 L (2,400 cubic in)
Warranty	2-year manufacturer warranty