

SCIENTIFIC INDUSTRIES SUGGESTED CATALOG COPY

Suggested Category: Incubators, Hybridization Oven/Incubator, Shakers, Rotators, Rockers

New!! Incubator-Genie™— Benchtop Shaking / Rotating Incubator

The multi-functional Incubator-Genie™ is designed to Incubate, Rotate, Rock, and Shake all in ONE compact unit.

- Precise temperature control and uniformity – range from 28°C to 75°C.
- Magnetic stainless steel platform may be easily set to rock or rotate almost any container.
- Containers can be easily mounted on the magnetic platform with magnetic clips or strips.
- Ideal for Hybridization.
- Chamber air is gently and continuously re-circulated ensuring temperature consistency and uniformity.
- A variety of accessories are available to provide maximum application versatility.

ORDERING INFORMATION

Catalog No.	Description
SI-1400	Incubator-Genie, 120V
SI-1401	Incubator-Genie, 230V, No Plug
0A-1200-001	Non-Skid Tray
SI-1120	6 Clip Plates for 12 each 10-13mm tubes
SI-1121	6 Clip Plates for 6 each 15-17mm tubes
SI-1122	6 Clip Plates for 3 each 28-30mm tubes
SI-1123	1 Clip Plate for 12 each 10-13mm tubes
SI-1124	1 Clip Plate for 6 each 15-17mm tubes
SI-1125	1 Clip Plate for 3 each 28-30mm tubes
SI-1126	Magnetic Bag Mounting Strips (pack of 16)
SI-1127	Expansion Kit (two stainless steel trays, 24 magnetic bag mounting strips)
SI-1130	2 Clip Plates for Hybridization Tubes
SI-1131	Stackable Shelf
SI-1134	2 Universal Clip Plates
SI-1135	1 Magnetic Covered 100 Microtube Box
SI-1136	1 Clip Plate for 100ml Volumetric Flasks
SI-1250	3-D Orbital Shaker Attachment
SI-1275	Horizontal Orbital Shaker Attachment

SPECIFICATIONS

Supplied with	Magnetic Platform, 1 Clip Plate each for 10-13mm, 15-17mm, and 28-30mm Tubes, 4 Magnetic Bag Mounting Strips, Non-Skid Tray, and Shelf.
Rotating	3-35 RPM
Rocking	6-70 cycles per minute

TEMPERATURE AND SETTINGS

Temperature Range	28°C - 75°C
Accuracy	±0.2°C
Uniformity	±0.5°C
Resolution	1.0°C
Stability	±0.1°C

DIMENSIONS

Chamber	254 x 362 x 260mm (D x W x H) (10 x 14.25 x 10.25in)
Overall	400 x 560 x 370mm (D x W x H) (15.75 x 22 x 14.5in)
Platform Capacity	4.5 Kg (10 lb)
Chamber Volume	24 L (0.86 ft³)
Weight	19.5 Kg (43 lb)
Warranty	2-year manufacturer warranty