

Suggested Category: **Stirrers**

## **MEGAMAG™ GENIE LR – Low Speed Single Position Large Volume Magnetic Stirrer**

Large volume magnetic stirrer designed for low speed large capacity applications.

- Stirring speed of 50-300 RPM.
- Powerful magnet and motor offer exceptional magnetic force capable of stirring up to 38L.
- Digital feedback (true speed) adjustable in 10 RPM increments.
- LED display for easier speed and timed readings.
- Rugged ABS/Polycarbonate housing and the famous Vortex-Genie heavy-duty motor provide years of reliable service.
- Spill-resistant housing channels fluids away from the internal components.
- Can be used in cold rooms and incubators.
- Non-Skid feet.

### **ORDERING INFORMATION**

SI-3236L	Digital MegaMag Genie LR, 120V
SI-3246L	Digital MegaMag Genie LR, 230V, No Plug

### **SPECIFICATIONS**

Stirring Speed	50-300 RPM
Time	1-99 min or continuous
Maximum Stirring Positions	1
Maximum Stirring Capacity	38 L
Stir Bar Maximum Length	7.6 cm (3in)
Material	ABS/Polycarbonate
Top Surface Dimensions (D x W x H)	305 x 305mm 12 x 12in
Overall Dimensions (D x W x H)	310 x 305 x 130mm 12.25 x 12 x 5.1in
Weight	3 Kg (6.6 lb)