

Suggested Category: Stirrers

MEGAMAG™ GENIE –Single Position Large Volume Magnetic Stirrer

Large volume magnetic stirrer designed for large capacity applications.

- Powerful magnet and motor offer exceptional magnetic force capable of stirring up to 38L.
- Both analog control and digital feedback (true speed) versions are available for MegaMag Genie.
- True Speed readout (RPM) with feedback control on digital version of the MegaMag genie is adjustable in 10 RPM increments.
- LED display is offered on digital version for easier speed and timed readings.
- Rugged ABS/Polycarbonate housing and the famous Vortex-Genie heavy-duty motor provide years of reliable service.
- Spill-resistant housing channels fluids away from the internal components.
- Can be used in cold rooms and incubators.
- Non-Skid feet.

ORDERING INFORMATION

Catalog No.	Description
SI-2236	MegaMag Genie 120V
SI-2246	MegaMag Genie 230V, No Plug
SI-3236	Digital MegaMag Genie 120V
SI-3246	Digital MegaMag Genie 230V, No Plug

SPECIFICATIONS

Stirring Speed	250-1300 RPM
Time	1-99 min or continuous (digital only)
Maximum Stirring Positions	1
Maximum Stirring Capacity	38 L
Stir Bar Maximum Length	7.6 cm (3in)
Material	ABS/Polycarbonate Nonskid Top Pad (1)
Top Surface Dimensions (D x W x H)	305 x 305mm 12 x 12in
Overall Dimensions (D x W x H)	310 x 305 x 130mm 12.25 x 12 x 5.1in
Weight	3 Kg (6.6 lb)