

# SCIENTIFIC INDUSTRIES SUGGESTED CATALOG COPY

Suggested Category: Mixers, Shakers, Vortexers

## **MICROPLATE GENIE® Single Position Microplate Mixer**

The MicroPlate Genie provides high speed action for increased mixing efficiency in every well.

- Aggressive mixing speed from 750-3200 RPM for complete and uniform mixing in any microplate format.
- Small vortexing orbit of 1.0mm for thorough mixing regardless of sample viscosity.
- Speed and orbit combine to offer true vortexing action in each well of a microplate.
- Accepts most microplate types within the recommendations of the Society for Biomolecular Screening (SBS), even in 384-well formats.
- Built-in timer operates from 1 – 60 minutes or continuous for hands-free use.

### **ORDERING INFORMATION**

<b>Catalog No.</b>	<b>Description</b>
SI-0400	MicroPlate Genie, 120V
SI-0401	MicroPlate Genie, 230V, No Plug

### **SPECIFICATIONS**

Speed	750-3200 RPM (750 – 2700 RPM for 230V models)
Time	1 – 60 minutes or continuous
Vortexing Orbit	1.0mm diameter
Maximum Capacity	1 microplate, any type with standard SBS footprint
Dimensions (D x W x H)	150 x 180 x 165mm
	5.75 x 7.0 x 6.5in
Weight	2.75 Kg (6.1 lb)