

# SCIENTIFIC INDUSTRIES SUGGESTED CATALOG COPY

## **NEW Mini-100 Orbital-Genie™— Benchtop Orbital Shaker**

Slow-shaking orbital shaker with a 11.4"x11.4" (290x290mm) non-slip, rubber-coated platform with spill-proof edges is ideal for gentle shaking applications, with no additional accessories necessary. Variable shaking speed up to 100RPM.

Double or even triple capacity with optional stacking tray sold separately.

### **ORDERING INFORMATION**

<b>Catalog No.</b>	<b>Description</b>
SI-M100	Orbital-Genie Shaker, 120V
SI-M101	Orbital-Genie Shaker, 230V No Plug
SI-M150	Stacking Tray for Mini-100

### **SPECIFICATIONS**

Speed	10 – 100 RPM
Orbit	19mm
Platform Dimensions (D x W)	290 x 290mm (11.4 x 11.4in)
Overall Dimensions (D x W)	350 x 350mm (14 x 14in)
Weight	3.5 Kg (7.8 lbs.)
Warranty	2-year manufacturer warranty