

# SCIENTIFIC INDUSTRIES SUGGESTED CATALOG COPY

## NEW Mini-300 Orbital-Genie™— Benchtop Orbital Shakers

Powerful standard size orbital shakers with a 12"x12" (305x305mm) platform, available in three different configurations. Variable shaking speed up to 300RPM.

- Mini-300 Orbital-Genie with Universal Ratcheting Clamps – provides the flexibility for using a mix of different size and shape vessels on a single platform. Unit is supplied with three ratcheting cushioned clamps. Extra ratcheting clamps for more than two rows are optional.
- Mini-300 Orbital-Genie with Adhering Mat provides a flat platform with adhering rubber mat perfect for use with open containers. Patented material is naturally adhering. Rinse clean to renew adhering properties.
- Mini-300 Orbital-Genie for Standard Flask Clamps is equipped with a platform for accepting various size flasks or beakers. Clamps sold separately.

Double or even triple capacity with optional stacking tray sold separately.

### ORDERING INFORMATION

| Catalog No. | Description   |
|-------------|---|
| SI-M1500    | Mini-300 Orbital-Genie Shaker, Ratcheting Clamps, 120V              |
| SI-M1600    | Mini-300 Orbital-Genie Shaker, Standard Flask Clamps, 120V          |
| SI-M1700    | Mini-300 Orbital-Genie Shaker, Adhering Mat, 120V                   |
| SI-M1501    | Mini-300 Orbital-Genie Shaker, Ratcheting Clamps, 230V, No Plug     |
| SI-M1601    | Mini-300 Orbital-Genie Shaker, Standard Flask Clamps, 230V, No Plug |
| SI-M1701    | Mini-300 Orbital-Genie Shaker, Adhering Mat, 230V, No Plug          |
| SI-1610     | 100 mL Flask Clamp (Platform accepts 16 of this size)               |
| SI-1611     | 125/150 mL Flask Clamp (Platform accepts 16 of this size)           |
| SI-1612     | 250 mL Flask Clamp (Platform accepts 16 of this size)               |
| SI-1613     | 500 mL Flask Clamp (Platform accepts 9 of this size)                |
| SI-1614     | 1 L Flask Clamp (Platform accepts 5 of this size)                   |
| SI-1615     | 2 L Flask Clamp (Platform accepts 3 of this size)                   |
| SI-1616     | Double-Sided Adhesive Tape Roll, 60in. x 1in. (1.5m x25.4mm)        |
| SI-M1617    | Adhering Mat, Blue, 11.8in. x 11.8in. (300mm x 300mm)               |
| SI-M1510    | Universal Ratcheting Clamp for Mini-300                             |

### SPECIFICATIONS

|                                |                              |
|--------------------------------|------------------------------|
| Speed                          | 10 – 300 RPM                 |
| Orbit                          | 19mm                         |
| Platform Dimensions<br>(D x W) | 305 x 305mm (12 x12in)       |
| Overall Dimensions<br>(D x W)  | 350 x 350mm (14 x 14in)      |
| Weight                         | 10.25 Kg (22.6 lbs.)         |
| Warranty                       | 2-year manufacturer warranty |