

SCIENTIFIC INDUSTRIES SUGGESTED CATALOG COPY

Suggested Category: **Mixers, Shakers, Vortexers**

NEW!! MULTI VORTEX-GENIE® — Large Capacity Multi Vessel Vortex Mixer

Designed to hold and provide true vortexing in almost any container up to 250ml with same reliability found in the original Vortex-Genie 2.

- Holds multiple test tubes, vials, or cylinders of various sizes (racks included) or up to six up to 250ml vessels or combination of both.
- Digital control and display of speed and time for reproducibility and accuracy, with true speed display.
- Gentle ramp up to speed, so open containers will not spill.
- Use for gentle, orbital-type shaking or more vigorous vortexing.
- Durable and heavy zinc die casting with specially designed feet prevents “walking”.
- May be used in cold rooms and incubators.

ORDERING INFORMATION

Catalog No.	Description
SI-M236	Multi Vortex-Genie®, 120V
SI-M246	Multi Vortex-Genie®, 230V, No Plug

SPECIFICATIONS

Supplied with	Two 12-Pack Multi-Tube Racks holding up to 24 tubes of any kind or material up to 29mm. Dia. (50ml centrifuge tube)
Speed	250-1000 rpm
Orbit	4mm dia (same as Vortex-Genie 2)
Timer	1-99 minutes or continuous
Dimensions	310 x 270 x 195 mm (D x W x H) (12.2 x 10.6 x 7.7 in)
Weight	11.3 Kg (25 lbs)
Warranty	2-year manufacturer warranty