

Suggested Category: Stirrers

QUADMAG™ GENIE- Four-Position Magnetic Stirrer

Powerful, general-purpose 4-position magnetic stirrer, designed for general laboratory applications with easy-to-use controls that allow users to easily change the desired speed.

- Powerful magnet and motor offers exceptional magnetic force capable of stirring up to 2 liters for each position.
- Powerful “rare earth” magnets automatically pull stir bars into position and stay coupled even under the most demanding stirring conditions
- Rugged ABS/Polycarbonate housing and the famous Vortex-Genie heavy-duty motor provide years of reliable service
- Spill-resistant housing channels fluids away from the internal components.
- Can be used in cold rooms and incubators.
- Non-skid feet.

Ordering Information

Catalog No.

SI-4236

SI-4246

Description

QuadMag Genie 120V

QuadMag Genie 230V, No Plug

Specifications

Stirring Speed	250-1300 RPM
Time	N/A
Maximum Stirring Positions	4
Maximum Stirring Capacity	2 L each position
Stir Bar Maximum Length	2in
Material	ABS/Polycarbonate Nonskid Top Pads (4)
Top Surface Dimensions (D x W x H)	305 x 305mm 12 x 12in
Overall Dimensions (D x W x H)	310 x 305 x 130mm 12.25 x 12 x 5.1in
Weight	3 Kg (6.6 lb)