Cell Disruptor



- Rapid hands-free disruption of up to 12 microtubes (1.5ml or 2.0ml) with NEW Universal Microtube Holder
- Timer for speed up to 15 minutes for Analog model, up to 99 minutes for Digital model or continuous for both models
- Suitable for use in cold rooms or incubators
- Less costly than ultrasonic units or homogenizers
- Use with Disruptor Beads™, made from spherical lead free soda lime glass, available in two sizes
- Disruptor Beads are used with any Disruptor Genie or TurboMix[™] attachment for Vortex-Genie[®] 2 family

Digital model is ideal for accurate, reproducible, repeatable results. Features adjustable speed and speed alarm.

Disruptor Genie and Disruptor Beads applications:

- Disruption of bacterial cultures for nucleic acid isolation and downstream processing
- Grinding plant tissue for DNA extractions
- Preparation of cell lysates from yeast or fungi
- Rapid resuspension of pellets
- Dissolution of chemical substances

	Ordering Information
Catalog No.	Description
SI-D238	Disruptor Genie, 120V
SI-D248	Disruptor Genie, 230V, No Plug
SI-D258	Disruptor Genie, 230V, Euro Plug
SI-D268	Disruptor Genie, 230V, British Plug
SI-D278	Disruptor Genie, 230V, Swiss Plug
SI-D298	Disruptor Genie, 240V, Australian Plug
SI-D288	Disruptor Genie, 100V
SI-DD38	Digital Disruptor Genie, 120V
SI-DD48	Digital Disruptor Genie, 230V, No Plug
SI-DD58	Digital Disruptor Genie, 230V, Euro Plug
SI-DD68	Digital Disruptor Genie, 230V, British Plug
SI-DD78	Digital Disruptor Genie, 230V, Swiss Plug
SI-DD98	Digital Disruptor Genie, 240V, Australian Plug
SI-DD88	Digital Disruptor Genie, 100V
SI-BG01	Disruptor Beads, 0.1mm, 375g (8 fl. oz.)
SI-BG05	Disruptor Beads, 0.5mm, 375g (8 fl. oz.)

Specifications

Weight 4.3 Kg (9.5 lb) **Base Dimensions** 165 x 122 x 190mm $(D \times W \times H)$ 6.5 x 4.8 x 7.5in Time range (Analog)

0 - 15 minutes or continuous Time range (Digital) 0 - 99 minutes or continuous Capacity Up to 12 - 1.5ml or 2.0ml microtubes

(additional holders sold separately, tubes not supplied)

3000 RPM (2850 for 50Hz models) Speed (Analog) Speed (Digital) 1000-3000 RPM (2850 for 50Hz models)

Speed Alarm (Digital) On/Off, Adjustable

> 0.1mm diameter beads for bacteria 0.5mm diameter beads for yeast

375 grams (8 fl. oz.)



(Analog Only)

ETL Listed Conforms to UL Std. 3101-1 ETL Listed Conforms to CA State ETL Listed Certified to CAN/CSA C22.2







Disruptor Beads Bead Sizes

Weight

Disruptor Genie



Home of the Multi-Task Genie™ Family

80 Orville Drive, Suite 102 Bohemia, NY 11716 U.S.A.

Phone: (631) 567-4700 • Fax: (631) 567-5896 Toll Free: 888-850-6208

Customer Service: custsvc@scientificindustries.com