

## Another revolutionary product from the makers of the Vortex-Genie®!

- O Holds a wide variety of tubes, flasks, microplates, bags, and odd-shaped vessels at any angle.
- O Available in two INTERCHANGEABLE platform configurations.
- O Standard platform is designed to turn 90 degrees allowing different mixing motions.
- MAGNETIC platform features unique spring clip plate system. (clip plates sold separately)
- O Designed for cold room and incubator use.
- Uniquely designed for a wide selection of motions adjustable from gentle rolling to forceful end over end mixing.
- O Continuous-duty quality motor provides smooth, quiet operation.

	Ordering Information	Specifications	
Catalog No. SI-2100 SI-2200 SI-2101 SI-2201 SI-2102 SI-2202 SI-2103	Description Rotator Genie, Standard Platform, 120V Rotator Genie, Magnetic Platform, 120V Rotator Genie, Standard Platform, 230V, No Plug Rotator Genie, Magnetic Platform, 230V, No Plug Rotator Genie, Standard Platform, 230V, Euro Plug Rotator Genie, Magnetic Platform, 230V, Euro Plug Rotator Genie, Standard Platform, 230V, British Plug	Supplied  Speed Overall Dimensions (W x D x H)  Platform Dimensions	Standard Platform with Retaining Bands or Magnetic Platform 5-35 RPM 400 x 175 (290 platform horizontal) x 400mm (15.8 x 7.0 (11.5 platform horizontal) x 15.8 in)
SI-2203 SI-2104 SI-2204 SI-2106 SI-2206 SI-2105 SI-2205	Rotator Genie, Magnetic Platform, 230V, British Plug Rotator Genie, Standard Platform, 230V, Swiss Plug Rotator Genie, Magnetic Platform, 230V, Swiss Plug Rotator Genie, Standard Platform, 240V, Aust. Plug Rotator Genie, Magnetic Platform, 240V, Aust. Plug Rotator Genie, Standard Platform, 100V Rotator Genie, Magnetic Platform, 100V	Standard Platform Magnetic Platform Weight Environmental	260 x 260mm (10.25 x 10.25 in) 290 x 315mm (11.5 x 12.5 in) 5.2 kgs (11.5 lbs) Installation Category 2, 0° - 38°C (32° - 100°F), 95% max. relative humidity









80 Orville Drive, Suite 102 Bohemia, NY 11716 U.S.A. scientificindustries.com

## **Rotator Genie™ Accessories** Fig. 2 Fig. 1 Fig. 12 Fig. 14 Fig. 11 Fig. 13 Fig. 9 Fig. 4, 7 Fig. 3, 6 Fig. 5, 8 Fig. 10

Fig. No.	Catalog No.	Description	
1	SI-2110	Standard Platform with Bands	
2	SI-2210	Magnetic Platform	
3	SI-1123(*)	1 Clip Plate for 12 each 10-13mm tubes (magnetic platform capacity: 12 clip plates)	
4	SI-1124(*)	1 Clip Plate for 6 each 15-17mm tubes (magnetic platform capacity: 12 clip plates)	
5	SI-1125(*)	1 Clip Plate for 3 each 28-30mm tubes (magnetic platform capacity: 12 clip plates)	
6	SI-1120(*)	6 Clip Plates for 12 each 10-13mm tubes (magnetic platform capacity: 12 clip plates)	
7	SI-1121(*)	6 Clip Plates for 6 each 15-17mm tubes (magnetic platform capacity: 12 clip plates)	
8	SI-1122(*)	6 Clip Plates for 3 each 28-30mm tubes (magnetic platform capacity: 12 clip plates)	
9	SI-1126(*)	Bag Mounting Strips (pack of 16)	
10	SI-1130(*)	2 Clip Plates for Hybridization Tubes	
11	SI-1134(*)	2 Universal Clip Plates	
12	SI-1135(*)	1 Magnetic Covered 100 Microtube Box	
13	SI-1136(*)	1 Clip Plate for 100ml Volumetric Flasks	
14	SI-1140(*)	1 Magnetic 50 Microtube Box	

<sup>(\*) -</sup> Requires magnetic platform, (Cat. No. SI-2210), which is supplied with all SI-22XX models or can be purchased separately.



Scientific Industries, Inc. Home of the Multi-Task Genie™ Family

80 Orville Drive, Suite 102 Bohemia, NY 11716 U.S.A. scientificindustries.com