
SITM NEWS

WHAT'S NEW AT...

Scientific Industries, Inc.
70 Orville Drive
Bohemia, New York 11716

FOR IMMEDIATE RELEASE (Bohemia, NY, February 9, 2004)

Scientific Industries Introduces New Glass Bead Media To Complement Existing Products

Scientific Industries Inc. (OTCBB-SCND) announces the availability of Disruptor BeadsTM borosilicate glass bead media. These spherical beads are commonly used by researchers for mechanical disruption of many sample types especially bacteria and yeast cells. The use of Disruptor Beads, along with Scientific Industries' Disruptor Genie[®] and TurboMixTM attachment for the Vortex-Genie[®] 2 or Vortex-Genie[®] 2T, may be used for rapid and effective sample disruption in up to 12 microtubes. The Disruptor Beads are available in two sizes (0.5 mm diameter and 0.25 mm diameter) and may be purchased from Scientific Industries at www.scientificindustries.com.

About Scientific Industries

Scientific Industries manufactures and markets laboratory equipment, including the world-renowned Vortex-Genie[®] 2 Mixer. Scientific Industries' products are used by research laboratories, clinics, pharmaceutical manufacturers, medical device manufacturers, and other industries.

“Statements made in this press release that relate to future events, performance or financial results of the Company are forward-looking statements which involve uncertainties that could cause actual events, performance or results to materially differ. The Company undertakes no obligation to update any of these statements. Readers are cautioned not to place undue reliance on these forward-looking statements, which speak only as to the date hereof. Accordingly, any forward-looking statement should be read in conjunction with the additional information about risks and uncertainties set forth in the Company's Securities and Exchange Commission reports, including our annual report on Form 10-KSB.”

###



Phone: 888-850-6208
Fax: 631-567-5896
info@scientificindustries.com
www.scientificindustries.com