



WHAT'S NEW AT...

Scientific Industries, Inc.

Airport International Plaza
70 Orville Drive

Bohemia, New York 11716

FOR IMMEDIATE RELEASE (Bohemia, NY, January 3, 2000)

Scientific Industries, Inc. Announced Today an Innovative New Product

ENVIRO-GENIE™

Your Own Personal Cold Room / Incubator

Scientific Industries, Inc. (OTCBB:SCND). The Company introduced the Enviro-Genie™, a unique **environmentally controlled benchtop 5-in-1 refrigerator / incubator / rotator / rocker / stirrer** designed to provide the security of precise temperature control ranging from 4°C to 75°C for demanding incubation/refrigeration applications without taking up a great deal of space. The Enviro-Genie™ is equipped with two revolutionary magnetic stirrers. The magnetic stirrers feature a speed range of 1-1000 RPM and include stir-reverse. Chamber air is gently and continuously circulating at a rate that ensures the temperature consistency and uniformity of all test samples. A shatter resistant viewing window allows for inspection of samples without affecting the chamber environment.

The Enviro-Genie's rotating / rocking platform is magnetized to accommodate a non-skid metal rocker platform or various sized autoclaveable spring clip mounting plates. It performs membrane-based assays / washings and hybridizations more efficiently. The magnetic platform can be easily removed and up to two stackable wire racks can be placed inside the chamber (.86 cu ft.) to be used as a standard refrigerator/incubator for staining trays, plastic bags, tissue culture plates, etc.

The Enviro-Genie microprocessor controls easily achieve desired rotating, rocking and stirring speed; temperature, and time; all of which are displayed on large LEDs for accurate reproducible conditions. Timer displays elapsed time or remaining time until programmed off - up to 96 hours. Visible and audible alarms for both temperature and time.

###