

The logo for SI NEWS, featuring the letters 'SI' in a stylized font followed by the word 'NEWS' in a larger, bold, sans-serif font.

***The Latest Information...***



**Scientific Industries, Inc.**  
70 Orville Drive  
Bohemia, New York 11716

FOR IMMEDIATE RELEASE (Bohemia, NY, November 24, 2008)

### **Scientific Industries Launches a New Product**

Scientific Industries, Inc. (OTCBB-SCND), November 24, 2008 – Scientific Industries Inc., announced today its market launch of the Incubator-Genie<sup>®</sup>, a new addition to its line of shakers and vortex mixers, including the world-renowned Vortex-Genie. The innovative and patented design of the Incubator-Genie<sup>®</sup>, a benchtop shaking/rotating incubator, incorporates a magnetic stainless steel platform which can be set for use as a rotator or rocker. The unit can be used with a variety of containers which are held on the platform using various magnetic clips and strips.

The new product, which is ideal for a variety of applications in research laboratories, will be offered world-wide through the Company's extensive network of laboratory equipment dealers.

#### **About Scientific Industries**

Scientific Industries designs, manufactures, and markets a variety of laboratory equipment, including the world-renowned Vortex-Genie<sup>®</sup> 2 Mixer, and since November 2006 produces and sells customized catalyst research instruments. Scientific Industries' products are generally used for research purposes in laboratories of universities, hospitals, pharmaceutical companies, chemical companies, and medical device manufacturers.

###



Phone: 888-850-6208

Fax: 631-567-5896

[info@scientificindustries.com](mailto:info@scientificindustries.com)

[www.scientificindustries.com](http://www.scientificindustries.com)